Cirona[™] 6200

Deep Vein Thrombosis Prevention Therapy System

Product Overview

The Cirona[™] 6200 is a compression therapy pump used for patients who may be at-risk of developing a painful and potentially life threatening blood clot known as deep vein thrombosis or DVT. The device comes with a pair of soft compression garments what wrap around the calf, calf-thigh, or foot on each leg or foot. Air compression inflates the garments during a 60 second automatically timed cycle (an approximate 12-second inflation period followed by a 48-second relaxation period) to mimic the body's natural walking functionality, to enhance circulation and prevent early blood clothing.

Features and Benefits:



Convenient comes with long-lasting rechargeable battery, bed hanger attaches easily in hospital or home care setting.



Effective: Enhances circulation to provide therapy that may prevent clotting.



Comfortable: Disposable compression garments mode of soft, breathable material.

Device Specifications:

Weight	5.1 lbs.	
Dimensions	4" (W) x 7" (L) x 11" (H)	
Inflation Cycle	12 seconds	
Rest Period	48 seconds	
Cycle Time	60 seconds	
Pressure	40mmHg	
Electrical	100 - 240 VAC, 50x60 Hz	
Fuse Rated	250 VAC, 1.0AMP, SLO-BLO	
Operation Mode	Continuous	
Battery :	Lithium Battery, 10.8V, 4400mAh (rechargeable)	

Marketed in India By:







Cirona[™] Sleeves

DeepVein Thrombosis Prevention Therapy System

Model	Application	Circumference	Patient Sizing
D-110	Calf	Up to 17"	Standard
D-120	Calf	Up to 23"	Large
D-125	Calf	Up to 32"	Extra Large
D-135	Calf-Thigh	Up to 22"	Small
D-130	Calf-Thigh	Up to 28"	Standard
D-140	Calf-Thigh	Up to 36"	Large
D-145	Calf-Thigh	Up to 42"	Extra Large
F-15	Foot	Up to 13"	Small
F-10	Foot	Up to 16"	Standard
Extension Tubing Set	6' tube extender to connect the garments to the device		

The Cirona Sleeves are compatible with the Cirona 6100 and Cirona 6200 Dee Vein Thrombosis Prevention Therapy Pumps. The 6300, our portable model of this device, comes with the sleeve attached.















