

 **ACCUREX**
— Since 1984 —

Xpress Alc

THE POWER OF HbA1c IN YOUR HAND

US FDA



HANDHELD HbA1c METER

BATTERY OPERATED

NO PIPETTING, NO LYSING

LAB LIKE ACCURACY





Unique

- Battery Operated
- Truly Portable
- Voice prompt for User Guidance



Accurate

- Accurate result with cv <3%
- Lab like accuracy with 95%
- No interference of HB Variant



Convenient

- No pipeting
- No lysing
- Reagent: Suitable to store at room temperature (upto 30° C)

“Product with new Technology, patented in 13 countries including India.”

Specification	Technical Parameters
Testing Principle	Boronate Affinity Chromatography
Testing Parameter	Glycohemoglobin (HbA1c)
Measuring Range	4.0% - 14.0%
Precision	CV <3% (HbA1c: 4.0% - 6.5%)
Blood Sample	Finger Prick or Venous Blood (EDTA Anticoagulant)
Blood Volume	About 3ul
Testing Time	About 5 minutes
Data Unit	Set in Advance the Data Unit • NGSP%; IFCC mmol/mol; eAG* mmol/l
Voice Prompt	Voice Prompt to Guide the user
Data Storage	1000 Testing Data
Data Port	Mini USB data interface can be connected with HIS/LIS; System Thermal Printer
Power Required	AAA Battery x 4
Analyser Dimension	61.6 x 122.9 x 24.5 mm
Screen Size	47 x 32 mm

* eAg conversion factor
mmol/ l X 18.0182 = mg/ dl



ACCUREX BIOMEDICAL PVT. LTD.

Head Office - Mumbai. Tel. :91(022) 6744 6744

Customer Support no.:1800 - 209 8456

E-mail: info@accurex.net Website: www.accurex.org

Plant : G-54, MIDC Tarapur, Boisar, Thane - 401506. India